

REFERRED_DB

PostgreSQL Database Schema

Source	Schema-only pg_dump provided by the user
Database engine	PostgreSQL 16.14
Schema	public
Tables	75
Foreign-key relationships	53

This document describes the database structure only. The supplied export contains definitions, constraints, indexes, and relationships, but no table data.

1. Database Overview

Source of truth. This document was generated from the PostgreSQL schema-only dump, so table names, columns, data types, keys, constraints, and foreign-key references are taken from the supplied database definition.

Tables	Columns	Primary keys	Foreign keys	Indexes
75	939	75	53	93

1.1 Notation

- **PK:** Primary key
- **FK:** Foreign key
- **UQ:** Single-column unique constraint
- **NN:** NOT NULL
- **IDENTITY:** PostgreSQL generated identity column

1.2 Organizational grouping

The categories below are inferred from table names and Django naming conventions; PostgreSQL itself does not store these application-group labels in the dump.

Group	Count	Tables
Core referred/final workflow	17	Antibiotic_List, AntibioticEntry, AntibioticEntry_ab_breakpoints_id, Batch_Table, BreakpointsTable, Classification_Table, Emerging_Filter_Age, Emerging_Table, Final_Data, FinalAntibioticEntry, FinalAntibioticEntry_ab_breakpoints_id, Phenotype_Post, Phenotype_Pre, Recommendation_items, Referred_Data, SpecimenTypeTable, TATform
WGS and pipeline data	15	amrfinderplus, assembly-scan, BactScout, builtin_wgs_pipeline_setting, checkm2, custom_wgs_pipeline, custom_wgs_pipeline_field, custom_wgs_pipeline_record, custom_wgs_pipeline_upload_batch, gambit, gtdbtk, mlst, sample_information, wgs_app_gambitdisplayconfig, WGS_Project
Upload/staging tables	22	AmrfinderUpload, Antibiotic_upload, AssemblyUpload, BactScoutUpload, Breakpoint_upload, Checkm2Upload, DemogsData_upload, Emerging_upload, FinalAntibiotic_upload, GambitUpload, GtdbTkUpload, home_pheno_upload_pre, home_reco_item_upload, MlstUpload, Organism_upload, Pheno_upload_Post, RawAntibiotic_upload, Referred_upload, SampleInfoUpload, SiteCode_upload, Specimen_upload, TATStepConfigUpload
Application support/reference tables	11	home_arsstaff_details, home_fieldmapping, home_final_concordancedetail, home_final_concordancereport, home_nonworkingday, home_organism_list, home_recommendation, home_sitedata, home_tatstep, home_tatstepconfig, UserProfile
Django/authentication system tables	10	auth_group, auth_group_permissions, auth_permission, auth_user, auth_user_groups, auth_user_user_permissions, django_admin_log, django_content_type, django_migrations, django_session

3. Foreign-Key Relationships

The schema contains 53 foreign-key column relationships.

Source table	Column	Referenced table	Referenced column	Constraint
amrfinderplus	amrfinder_project_id	WGS_Project	id	amrfinderplus_amrfinder_project_id_ac8762d1_fk_WGS_Project_id
AntibioticEntry	ab_idNum_referred_id	Referred_Data	AccessionNo	AntibioticEntry_ab_idNum_referred_id_24c049c7_fk_Referred_
AntibioticEntry_ab_breakpoints_id	antibioticentry_id	AntibioticEntry	id	AntibioticEntry_ab_b_antibioticentry_id_a8615639_fk_Antibioti
AntibioticEntry_ab_breakpoints_id	breakpointstable_id	BreakpointsTable	id	AntibioticEntry_ab_b_breakpointstable_id_4ca37068_fk_Breakpoin
assembly-scan	assembly_project_id	WGS_Project	id	assembly-scan_assembly_project_id_c76c862a_fk_WGS_Project_id
auth_group_permissions	group_id	auth_group	id	auth_group_permissions_group_id_b120cbf9_fk_auth_group_id
auth_group_permissions	permission_id	auth_permission	id	auth_group_permission_permission_id_84c5c92e_fk_auth_perm
auth_permission	content_type_id	django_content_type	id	auth_permission_content_type_id_2f476e4b_fk_django_co
auth_user_groups	group_id	auth_group	id	auth_user_groups_group_id_97559544_fk_auth_group_id
auth_user_groups	user_id	auth_user	id	auth_user_groups_user_id_6a12ed8b_fk_auth_user_id
auth_user_user_permissions	permission_id	auth_permission	id	auth_user_user_permission_permission_id_1fbb5f2c_fk_auth_perm
auth_user_user_permissions	user_id	auth_user	id	auth_user_user_permissions_user_id_a95ead1b_fk_auth_user_id
BactScout	bactscout_project_id	WGS_Project	id	BactScout_bactscout_project_id_c5aa94f6_fk_WGS_Project_id
Batch_Table	created_by_id	auth_user	id	Batch_Table_created_by_id_53bb9736_fk_auth_user_id
BreakpointsTable	Abx_List_Whonet_Abx	Antibiotic_List	Whonet_Abx	BreakpointsTable_Abx_List_Whonet_Abx_2e9bc96e_fk_Antibioti
checkm2	checkm2_project_id	WGS_Project	id	checkm2_checkm2_project_id_b6750be1_fk_WGS_Project_id
Classification_Table	Class_idNumReferred_id	Final_Data	id	Classification_Table_Class_idNumReferred__3f17e0d1_fk_Final_Dat
custom_wgs_pipeline	created_by_id	auth_user	id	custom_wgs_pipeline_created_by_id_af63e8f2_fk_auth_user_id
custom_wgs_pipeline_field	pipeline_id	custom_wgs_pipeline	id	custom_wgs_pipeline_pipeline_id_360558a9_fk_custom_wg
custom_wgs_pipeline_record	matched_final_data_id	Final_Data	f_AccessionNo	custom_wgs_pipeline_matched_final_data_id_13244890_fk_Final_Dat
custom_wgs_pipeline_record	matched_raw_data_id	Referred_Data	AccessionNo	custom_wgs_pipeline_matched_raw_data_id_254a30fb_fk_Referred_
custom_wgs_pipeline_record	pipeline_id	custom_wgs_pipeline	id	custom_wgs_pipeline_pipeline_id_b1ddbc96_fk_custom_wg
custom_wgs_pipeline_record	upload_batch_id	custom_wgs_pipeline_upload_batch	id	custom_wgs_pipeline_upload_batch_id_b8262770_fk_custom_wg
custom_wgs_pipeline_record	uploaded_by_id	auth_user	id	custom_wgs_pipeline_uploaded_by_id_921db27c_fk_auth_us
custom_wgs_pipeline_upload_batch	pipeline_id	custom_wgs_pipeline	id	custom_wgs_pipeline_pipeline_id_ed554286_fk_custom_wg
custom_wgs_pipeline_upload_batch	uploaded_by_id	auth_user	id	custom_wgs_pipeline_uploaded_by_id_c4d361c7_fk_auth_us
django_admin_log	content_type_id	django_content_type	id	django_admin_log_content_type_id_c4bce8eb_fk_django_co

REFERRED_DB Database Schema

Source table	Column	Referenced table	Referenced column	Constraint
django_admin_log	user_id	auth_user	id	django_admin_log_user_id_c564eba6_fk_auth_user_id
Emerging_Table	eme_primary_key_id	Final_Data	id	Emerging_Table_eme_primary_key_id_e0f72bab_fk_Final_Data_id
Final_Data	f_Batch_id_id	Batch_Table	id	Final_Data_f_Batch_id_id_e5f4c5cf_fk_Batch_Table_id
Final_Data	f_Spec_Type_id	SpecimenTypeTable	id	Final_Data_f_Spec_Type_id_aa f41bd6_fk_SpecimenTypeTable_id
FinalAntibioticEntry	ab_idNum_f_referred_id	Final_Data	id	FinalAntibioticEntry_ab_idNum_f_referred_23352223_fk_Final_Dat
FinalAntibioticEntry_ab_breakpoints_id	breakpointstable_id	BreakpointsTable	id	FinalAntibioticEntry_breakpointstable_id_b2298850_fk_Breakpoi
FinalAntibioticEntry_ab_breakpoints_id	final_antibioticentry_id	FinalAntibioticEntry	id	FinalAntibioticEntry_final_antibioticentr_541d5d3e_fk_FinalAnti
gambit	gambit_project_id	WGS_Project	id	gambit_gambit_project_id_60b79a1a_fk_WGS_Project_id
gtdbtk	gtdbtk_project_id	WGS_Project	id	gtdbtk_gtdbtk_project_id_76e0ec8e_fk_WGS_Project_id
home_arsstaff_details	User_Account_id	auth_user	id	home_arsstaff_details_User_Account_id_f48243d2_fk_auth_us
home_fieldmapping	user_id	auth_user	id	home_fieldmapping_user_id_6111f8c2_fk_auth_user_id
home_final_concordancedetail	report_id	home_final_concordancereport	id	home_final_concordance_report_id_428d5fa5_fk_home_fina
home_final_concordancereport	batch_id	Batch_Table	id	home_final_concordance_batch_id_d6ffb20_fk_Batch_Tab
home_final_concordancereport	created_by_id	auth_user	id	home_final_concordance_created_by_id_d9a0cf3a_fk_auth_user
home_final_concordancereport	final_data_id	Final_Data	id	home_final_concordance_final_data_id_466017ec_fk_Final_Dat
home_tatstep	performed_by_id	home_arsstaff_details	id	home_tatstep_performed_by_id_598e377c_fk_home_arss
home_tatstep	step_config_id	home_tatstepconfig	id	home_tatstep_step_config_id_6d1898ea_fk_home_tatstepconfi
home_tatstep	tat_id	TATform	id	home_tatstep_tat_id_2d827def_fk_TATform_id
mlst	mlst_project_id	WGS_Project	id	mlst_mlst_project_id_24fb710b_fk_WGS_Project_id
Referred_Data	Batch_id_id	Batch_Table	id	Referred_Data_Batch_id_id_2c0dc7b7_fk_Batch_Table_id
Referred_Data	Spec_Type_id	SpecimenTypeTable	id	Referred_Data_Spec_Type_id_35f213b0_fk_SpecimenTypeTable_id
sample_information	sample_project_id	WGS_Project	id	sample_information_sample_project_id_7ad83fbd_fk_WGS_Pr
SpecimenTypeTable	Specimen_Code_Grp_id	SpecimenTypeTable	id	SpecimenTypeTable_Specimen_Code_Grp_id_22b1f812_fk_Sp
TATform	tat_Batch_Isolates_id	Batch_Table	id	TATform_tat_Batch_Isolates_id_464d7e8d_fk_Batch_Table_id
UserProfile	user_id	auth_user	id	UserProfile_user_id_d5588aad_fk_auth_user_id
WGS_Project	Ref_Accession_id	Final_Data	f_AccessionNo	WGS_Project_Ref_Accession_id_0dcf3ffa_fk_Final_Dat

4. Table Catalog

The following catalog lists every table in the supplied public schema and its column count.

#	Table	Columns	Foreign keys out
1	amrfinderplus	27	1
2	AmrfinderUpload	2	0
3	Antibiotic_List	17	0
4	Antibiotic_upload	2	0
5	AntibioticEntry	46	1
6	AntibioticEntry_ab_breakpoints_id	3	2
7	assembly-scan	28	1
8	AssemblyUpload	2	0
9	auth_group	2	0
10	auth_group_permissions	3	2
11	auth_permission	4	1
12	auth_user	11	0
13	auth_user_groups	3	2
14	auth_user_user_permissions	3	2
15	BactScout	64	1
16	BactScoutUpload	2	0
17	Batch_Table	28	1
18	Breakpoint_upload	2	0
19	BreakpointsTable	25	1
20	builtin_wgs_pipeline_setting	7	0
21	checkm2	18	1
22	Checkm2Upload	2	0
23	Classification_Table	13	1
24	custom_wgs_pipeline	17	1
25	custom_wgs_pipeline_field	13	1
26	custom_wgs_pipeline_record	14	5
27	custom_wgs_pipeline_upload_batch	14	2
28	DemogsData_upload	2	0
29	django_admin_log	8	2
30	django_content_type	3	0
31	django_migrations	4	0
32	django_session	3	0
33	Emerging_Filter_Age	2	0
34	Emerging_Table	21	1
35	Emerging_upload	2	0
36	Final_Data	94	2
37	FinalAntibiotic_upload	2	0
38	FinalAntibioticEntry	47	1
39	FinalAntibioticEntry_ab_breakpoints_id	3	2
40	gambit	15	1
41	GambitUpload	2	0
42	gtdbtk	24	1
43	GtdbTkUpload	2	0
44	home_arsstaff_details	11	1
45	home_fieldmapping	7	1
46	home_final_concordancedetail	12	1
47	home_final_concordancereport	20	3
48	home_nonworkingday	5	0
49	home_organism_list	16	0
50	home_pheno_upload_pre	2	0
51	home_reco_item_upload	2	0
52	home_recommendation	3	0
53	home_sitedata	19	0
54	home_tatstep	11	3
55	home_tatstepconfig	5	0
56	mlst	14	1
57	MlstUpload	2	0

REFERRED_DB Database Schema

#	Table	Columns	Foreign keys out
58	Organism_upload	2	0
59	Pheno_upload_Post	2	0
60	Phenotype_Post	2	0
61	Phenotype_Pre	2	0
62	RawAntibiotic_upload	2	0
63	Recommendation_items	3	0
64	Referred_Data	107	2
65	Referred_upload	2	0
66	sample_information	25	1
67	SampleInfoUpload	2	0
68	SiteCode_upload	2	0
69	Specimen_upload	2	0
70	SpecimenTypeTable	6	1
71	TATform	19	1
72	TATStepConfigUpload	2	0
73	UserProfile	3	1
74	wgs_app_gambitdisplayconfig	3	0
75	WGS_Project	18	1

5. Detailed Table Definitions

For each table, the column table records the PostgreSQL type, nullability, key status, and foreign-key target. Constraint and index summaries follow each definition.

5.1 amrfinderplus

27 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Amrfinder_Accession	character varying(255)	YES			
name	character varying(255)	YES			
protein_id	character varying(255)	YES			
contig_id	character varying(255)	YES			
start	character varying(255)	YES			
stop	character varying(255)	YES			
strand	character varying(255)	YES			
element_symbol	character varying(255)	YES			
element_name	character varying(255)	YES			
scope	character varying(255)	YES			
type_field	character varying(255)	YES			
subtype	character varying(255)	YES			
class_field	character varying(255)	YES			
subclass	character varying(255)	YES			
method	character varying(255)	YES			
target_length	character varying(255)	YES			
reference_sequence_length	character varying(255)	YES			
percent_coverage_of_reference	character varying(255)	YES			
percent_identity_to_reference	character varying(255)	YES			
alignment_length	character varying(255)	YES			
closest_reference_accession	character varying(255)	YES			
closest_reference_name	character varying(255)	YES			
hmm_accession	character varying(255)	YES			
hmm_description	character varying(255)	YES			
Date_uploaded_am	date	NO			
amrfinder_project_id	bigint	YES	FK		WGS_Project.id

Constraints: amrfinderplus_pkey: PRIMARY KEY (id) | amrfinderplus_amrfinder_project_id_ac8762d1_fk_WGS_Project_id: FOREIGN KEY (amrfinder_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: amrfinderplus_Amrfinder_Accession_cb68fd3b (btree: "Amrfinder_Accession") | amrfinderplus_Amrfinder_Accession_cb68fd3b_like (btree: "Amrfinder_Accession" varchar_pattern_ops) | amrfinderplus_amrfinder_project_id_ac8762d1 (btree: amrfinder_project_id)

5.2 AmrfinderUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Amrfinderfile	character varying(100)	YES			

Constraints: AmrfinderUpload_pkey: PRIMARY KEY (id)

5.3 Antibiotic_List

17 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Show	boolean	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
Retest	boolean	NO			
Show_Site	boolean	NO			
Show_Ars	boolean	NO			
Show_Value	boolean	NO			
Disk_Abx	boolean	NO			
Tier	character varying(10)	NO			
Test_Method	character varying(100)	NO			
Abx_code	character varying(100)	NO			
Whonet_Abx	character varying(100)	NO	UQ		
Antibiotic	character varying(100)	NO			
Guidelines	character varying(100)	NO			
Potency	character varying(100)	NO			
Class	character varying(100)	NO			
Subclass	character varying(100)	NO			
Date_Modified	date	YES			

Constraints: Antibiotic_List_Whonet_Abx_key: UNIQUE ("Whonet_Abx") | Antibiotic_List_pkey: PRIMARY KEY (id)

Indexes: Antibiotic_List_Whonet_Abx_083a1de0_like (btree: "Whonet_Abx" varchar_pattern_ops)

5.4 Antibiotic_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_uploadAbx	character varying(100)	YES			

Constraints: Antibiotic_upload_pkey: PRIMARY KEY (id)

5.5 AntibioticEntry

46 columns; 2 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
ab_Site_Org	character varying(100)	YES			
ab_Ret_Org	character varying(100)	YES			
ab_Disk_Abx	boolean	NO			
ab_Ret_Disk_Abx	boolean	NO			
ab_AccessionNo	character varying(100)	YES			
ab_RefNo	character varying(100)	YES			
ab_Antibiotic	character varying(100)	YES			
ab_Abx_code	character varying(100)	YES			
ab_Abx	character varying(100)	YES			
ab_Disk_value	integer	YES			
ab_Disk_RIS	character varying(4)	NO			
ab_Disk_enRIS	character varying(4)	NO			
ab_MIC_operand	character varying(4)	YES			
ab_MIC_value	numeric(10,5)	YES			
ab_MIC_RIS	character varying(4)	NO			
ab_MIC_enRIS	character varying(4)	NO			
ab_AlertMIC	boolean	NO			
ab_Alert_val	character varying(30)	YES			
ab_R_breakpoint	character varying(10)	YES			
ab_I_breakpoint	character varying(10)	YES			
ab_SDD_breakpoint	character varying(10)	YES			
ab_S_breakpoint	character varying(10)	YES			
ab_Org_Flag	boolean	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
ab_Abx_Flag	boolean	NO			
ab_Abx_Phenotype	character varying(100)	YES			
ab_Retest_Antibiotic	character varying(100)	YES			
ab_Retest_Abx_code	character varying(100)	YES			
ab_Retest_Abx	character varying(100)	YES			
ab_Retest_DiskValue	integer	YES			
ab_Retest_Disk_RIS	character varying(4)	NO			
ab_Retest_Disk_enRIS	character varying(4)	NO			
ab_Retest_MIC_operand	character varying(4)	YES			
ab_Retest_MICValue	numeric(10,5)	YES			
ab_Retest_MIC_RIS	character varying(4)	NO			
ab_Retest_MIC_enRIS	character varying(4)	NO			
ab_Retest_AlertMIC	boolean	NO			
ab_Retest_Alert_val	character varying(30)	YES			
ab_Ret_R_breakpoint	character varying(10)	YES			
ab_Ret_I_breakpoint	character varying(10)	YES			
ab_Ret_SDD_breakpoint	character varying(10)	YES			
ab_Ret_S_breakpoint	character varying(10)	YES			
ab_MICJoined	character varying(7)	YES			
ab_idNum_referred_id	character varying(255)	YES	FK		Referred_Data.AccessionNo
ab_Abx_Phenotype_Other	character varying(100)	YES			
ab_Date_uploaded_rd	date	NO			

Constraints: AntibioticEntry_pkey: PRIMARY KEY (id) | AntibioticEntry_ab_idNum_referred_id_24c049c7_fk_Referred_: FOREIGN KEY ("ab_idNum_referred_id") REFERENCES public."Referred_Data"("AccessionNo") DEFERRABLE INITIALLY DEFERRED

Indexes: AntibioticEntry_ab_idNum_referred_id_24c049c7 (btree: "ab_idNum_referred_id") | AntibioticEntry_ab_idNum_referred_id_24c049c7_like (btree: "ab_idNum_referred_id" varchar_pattern_ops)

5.6 AntibioticEntry_ab_breakpoints_id

3 columns; 4 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
antibioticentry_id	bigint	NO	FK		AntibioticEntry.id
breakpointstable_id	bigint	NO	FK		BreakpointsTable.id

Constraints: AntibioticEntry_ab_break_antibioticentry_id_break_33730a43_uniq: UNIQUE (antibioticentry_id, breakpointstable_id) | AntibioticEntry_ab_breakpoints_id_pkey: PRIMARY KEY (id) | AntibioticEntry_ab_b_antibioticentry_id_a8615639_fk_Antibioti: FOREIGN KEY (antibioticentry_id) REFERENCES public."AntibioticEntry"(id) DEFERRABLE INITIALLY DEFERRED | AntibioticEntry_ab_b_breakpointstable_id_4ca37068_fk_Breakpoin: FOREIGN KEY (breakpointstable_id) REFERENCES public."BreakpointsTable"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: AntibioticEntry_ab_breakpoints_id_antibioticentry_id_a8615639 (btree: antibioticentry_id) | AntibioticEntry_ab_breakpoints_id_breakpointstable_id_4ca37068 (btree: breakpointstable_id)

5.7 assembly-scan

28 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Assembly_Accession	character varying(255)	YES			
sample	character varying(255)	YES			
total_contig	character varying(255)	YES			
total_contig_length	character varying(255)	YES			
max_contig_length	character varying(255)	YES			
mean_contig_length	character varying(255)	YES			
median_contig_length	character varying(255)	YES			
min_contig_length	character varying(255)	YES			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
n50_contig_length	character varying(255)	YES			
l50_contig_count	character varying(255)	YES			
num_contig_non_acgtn	character varying(255)	YES			
contig_percent_a	character varying(255)	YES			
contig_percent_c	character varying(255)	YES			
contig_percent_g	character varying(255)	YES			
contig_percent_t	character varying(255)	YES			
contig_percent_n	character varying(255)	YES			
contig_non_acgtn	character varying(255)	YES			
contigs_greater_1m	character varying(255)	YES			
contigs_greater_100k	character varying(255)	YES			
contigs_greater_10k	character varying(255)	YES			
contigs_greater_1k	character varying(255)	YES			
percent_contigs_greater_1m	character varying(255)	YES			
percent_contigs_greater_100k	character varying(255)	YES			
percent_contigs_greater_10k	character varying(255)	YES			
percent_contigs_greater_1k	character varying(255)	YES			
Date_uploaded_as	date	NO			
assembly_project_id	bigint	YES	FK		WGS_Project.id

Constraints: assembly-scan_pkey: PRIMARY KEY (id) | assembly-scan_assembly_project_id_c76c862a_fk_WGS_Project_id: FOREIGN KEY (assembly_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: assembly-scan_Assembly_Accession_b68fcab9 (btree: "Assembly_Accession") | assembly-scan_Assembly_Accession_b68fcab9_like (btree: "Assembly_Accession" varchar_pattern_ops) | assembly-scan_assembly_project_id_c76c862a (btree: assembly_project_id)

5.8 AssemblyUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Assemblyfile	character varying(100)	YES			

Constraints: AssemblyUpload_pkey: PRIMARY KEY (id)

5.9 auth_group

2 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	integer	NO	PK	YES	
name	character varying(150)	NO	UQ		

Constraints: auth_group_name_key: UNIQUE (name) | auth_group_pkey: PRIMARY KEY (id)

Indexes: auth_group_name_a6ea08ec_like (btree: name varchar_pattern_ops)

5.10 auth_group_permissions

3 columns; 4 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
group_id	integer	NO	FK		auth_group.id
permission_id	integer	NO	FK		auth_permission.id

Constraints: auth_group_permissions_group_id_permission_id_0cd325b0_uniq: UNIQUE (group_id, permission_id) | auth_group_permissions_pkey: PRIMARY KEY (id) | auth_group_permission_id_84c5c92e_fk_auth_perm: FOREIGN KEY (permission_id) REFERENCES public.auth_permission(id) DEFERRABLE INITIALLY DEFERRED | auth_group_permissions_group_id_b120cbf9_fk_auth_group_id: FOREIGN KEY (group_id) REFERENCES public.auth_group(id) DEFERRABLE INITIALLY DEFERRED

Indexes: auth_group_permissions_group_id_b120cbf9 (btree: group_id) | auth_group_permissions_permission_id_84c5c92e (btree: permission_id)

5.11 auth_permission

4 columns; 3 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	integer	NO	PK	YES	
name	character varying(255)	NO			
content_type_id	integer	NO	FK		django_content_type.id
codename	character varying(100)	NO			

Constraints: auth_permission_content_type_id_codename_01ab375a_uniq: UNIQUE (content_type_id, codename) | auth_permission_pkey: PRIMARY KEY (id) | auth_permission_content_type_id_2f476e4b_fk_django_co: FOREIGN KEY (content_type_id) REFERENCES public.django_content_type(id) DEFERRABLE INITIALLY DEFERRED

Indexes: auth_permission_content_type_id_2f476e4b (btree: content_type_id)

5.12 auth_user

11 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	integer	NO	PK	YES	
password	character varying(128)	NO			
last_login	timestamp with time zone	YES			
is_superuser	boolean	NO			
username	character varying(150)	NO	UQ		
first_name	character varying(150)	NO			
last_name	character varying(150)	NO			
email	character varying(254)	NO			
is_staff	boolean	NO			
is_active	boolean	NO			
date_joined	timestamp with time zone	NO			

Constraints: auth_user_pkey: PRIMARY KEY (id) | auth_user_username_key: UNIQUE (username)

Indexes: auth_user_username_6821ab7c_like (btree: username varchar_pattern_ops)

5.13 auth_user_groups

3 columns; 4 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
user_id	integer	NO	FK		auth_user.id
group_id	integer	NO	FK		auth_group.id

Constraints: auth_user_groups_pkey: PRIMARY KEY (id) | auth_user_groups_user_id_group_id_94350c0c_uniq: UNIQUE (user_id, group_id) | auth_user_groups_group_id_97559544_fk_auth_group_id: FOREIGN KEY (group_id) REFERENCES public.auth_group(id) DEFERRABLE INITIALLY DEFERRED | auth_user_groups_user_id_6a12ed8b_fk_auth_user_id: FOREIGN KEY (user_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: auth_user_groups_group_id_97559544 (btree: group_id) | auth_user_groups_user_id_6a12ed8b (btree: user_id)

5.14 auth_user_user_permissions

3 columns; 4 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
user_id	integer	NO	FK		auth_user.id
permission_id	integer	NO	FK		auth_permission.id

Constraints: auth_user_user_permissions_pkey: PRIMARY KEY (id) | auth_user_user_permissions_user_id_permission_id_14a6b632_uniq: UNIQUE (user_id, permission_id) | auth_user_user_permi_permission_id_1fbb5f2c_fk_auth_perm: FOREIGN KEY (permission_id) REFERENCES public.auth_permission(id) DEFERRABLE INITIALLY DEFERRED | auth_user_user_permissions_user_id_a95ead1b_fk_auth_user_id: FOREIGN KEY (user_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: auth_user_user_permissions_permission_id_1fbb5f2c (btree: permission_id) | auth_user_user_permissions_user_id_a95ead1b (btree: user_id)

5.15 BactScout

64 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
BactScout_Accession	character varying(255)	YES			
name	character varying(255)	YES			
status	character varying(255)	YES			
completeness	double precision	YES			
contamination	double precision	YES			
completeness_model_used	character varying(255)	YES			
translation_table_used	integer	YES			
coding_density	double precision	YES			
contig_n50	integer	YES			
average_gene_length	double precision	YES			
genome_size	integer	YES			
checkm2_gc_content	double precision	YES			
total_coding_sequences	integer	YES			
total_contigs	integer	YES			
max_contig_length	integer	YES			
additional_notes	text	YES			
a_final_status	character varying(255)	YES			
adapter_detection_status	character varying(255)	YES			
contamination_status	character varying(255)	YES			
species_status	character varying(255)	YES			
coverage_status	character varying(255)	YES			
coverage_estimate_qualibact_status	character varying(255)	YES			
duplication_status	character varying(255)	YES			
gc_content_status	character varying(255)	YES			
mlst_status	character varying(255)	YES			
n_content_status	character varying(255)	YES			
read_length_status	character varying(255)	YES			
read_q30_status	character varying(255)	YES			
species	character varying(255)	YES			
species_abundance	text	YES			
species_coverage	text	YES			
species_message	text	YES			
contamination_message	text	YES			
coverage_estimate_sylph	double precision	YES			
coverage_estimate_sylph_message	text	YES			
coverage_estimate_qualibact	double precision	YES			
coverage_estimate_qualibact_message	text	YES			
duplication_rate	double precision	YES			
duplication_message	text	YES			
gc_content	double precision	YES			
gc_content_lower	integer	YES			
gc_content_upper	integer	YES			
gc_content_message	text	YES			
n_content_rate	double precision	YES			
n_content_message	text	YES			
mlst_st	character varying(255)	YES			
mlst_message	text	YES			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
read1_mean_length	integer	YES			
read2_mean_length	integer	YES			
read_length_message	text	YES			
read_q20_bases	double precision	YES			
read_q20_rate	double precision	YES			
read_q30_bases	double precision	YES			
read_q30_rate	double precision	YES			
read_q30_message	text	YES			
read_total_bases	double precision	YES			
read_total_reads	integer	YES			
adapter_detection_message	text	YES			
ref_genome	character varying(255)	YES			
genome_size_expected	integer	YES			
Date_uploaded_bs	date	NO			
bactscout_project_id	bigint	YES	FK		WGS_Project.id
genome_size_expected_status	text	YES			

Constraints: BactScout_pkey: PRIMARY KEY (id) | BactScout_bactscout_project_id_c5aa94f6_fk_WGS_Project_id: FOREIGN KEY (bactscout_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: BactScout_BactScout_Accession_d9bbabea (btree: "BactScout_Accession") | BactScout_BactScout_Accession_d9bbabea_like (btree: "BactScout_Accession" varchar_pattern_ops) | BactScout_bactscout_project_id_c5aa94f6 (btree: bactscout_project_id)

5.16 BactScoutUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
bactscoutfile	character varying(100)	YES			

Constraints: BactScoutUpload_pkey: PRIMARY KEY (id)

5.17 Batch_Table

28 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
bat_SiteCode	character varying(255)	NO			
bat_Seq_No	integer	YES			
bat_Site_Name	character varying(255)	NO			
bat_Site_NameGen	character varying(255)	NO			
bat_Batch_Name	character varying(255)	NO			
bat_Batch_Code	character varying(255)	NO			
bat_Date_of_Entry	timestamp with time zone	NO			
bat_RefNo	character varying	YES			
bat_BatchNo	character varying(255)	NO			
bat_Total_batch	character varying(100)	NO			
bat_AccessionNo	character varying(255)	NO			
bat_AccessionNoGen	character varying(100)	NO			
bat_Default_Year	date	YES			
bat_Referral_Date	date	YES			
bat_Status	character varying(100)	NO			
bat_Encoder	character varying(255)	NO			
bat_Enc_Lic	character varying(100)	NO			
bat_Checker	character varying(255)	NO			
bat_Chec_Lic	character varying(100)	NO			
bat_Verifier	character varying(255)	NO			
bat_Ver_Lic	character varying(100)	NO			

Column	PostgreSQL type	Null	Key	Identity	Reference
bat_LabManager	character varying(255)	NO			
bat_Lab_Lic	character varying(100)	NO			
bat_Head	character varying(255)	NO			
bat_Head_Lic	character varying(100)	NO			
bat_Date_Accomplished	date	YES			
created_by_id	integer	YES	FK		auth_user.id

Constraints: Batch_Table_pkey: PRIMARY KEY (id) | Batch_Table_created_by_id_53bb9736_fk_auth_user_id: FOREIGN KEY (created_by_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: Batch_Table_created_by_id_53bb9736 (btree: created_by_id)

5.18 Breakpoint_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_uploadBP	character varying(100)	YES			

Constraints: Breakpoint_upload_pkey: PRIMARY KEY (id)

5.19 BreakpointsTable

25 columns; 2 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Guidelines	character varying(100)	NO			
Year	character varying(100)	NO			
Org	character varying(100)	NO			
Spec_code	character varying(100)	YES			
Emerging_specimen	boolean	NO			
Test_Method	character varying(20)	NO			
Potency	character varying(20)	NO			
Abx_code	character varying(15)	NO			
Tier	character varying(10)	NO			
Emerging_Org_Flag	boolean	NO			
Emerging_Abx_Flag	boolean	NO			
Emerging_Pheno_Flag	character varying(100)	NO			
Antibiotic	character varying(100)	NO			
Whonet_Abx	character varying(100)	NO			
Disk_Abx	boolean	NO			
R_val	character varying(30)	NO			
I_val	character varying(30)	NO			
SDD_val	character varying(30)	NO			
S_val	character varying(30)	NO			
Alert_val	character varying(30)	NO			
Date_Modified	date	NO			
Abx_List_Whonet_Abx	character varying(100)	YES	FK		Antibiotic_List.Whonet_Abx
Emerging_Pheno_Flag_Other	character varying(100)	NO			
Org_Code_type	character varying(100)	NO			

Constraints: BreakpointsTable_pkey: PRIMARY KEY (id) | BreakpointsTable_Abx_List_Whonet_Abx_2e9bc96e_fk_Antibioti: FOREIGN KEY ("Abx_List_Whonet_Abx") REFERENCES public."Antibiotic_List"("Whonet_Abx") DEFERRABLE INITIALLY DEFERRED

Indexes: BreakpointsTable_Abx_List_Whonet_Abx_2e9bc96e (btree: "Abx_List_Whonet_Abx") | BreakpointsTable_Abx_List_Whonet_Abx_2e9bc96e_like (btree: "Abx_List_Whonet_Abx" varchar_pattern_ops)

5.20 builtin_wgs_pipeline_setting

7 columns; 3 named constraints; 1 indexes.

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
pipeline_key	character varying(50)	NO	UQ		
display_name	character varying(100)	NO			
show_in_upload_center	boolean	NO			
show_in_overview	boolean	NO			
sort_order	integer	NO			
updated_at	timestamp with time zone	NO			

Constraints: builtin_wgs_pipeline_setting_sort_order_check: CHECK ((sort_order >= 0)) | builtin_wgs_pipeline_setting_pipeline_key_key: UNIQUE (pipeline_key) | builtin_wgs_pipeline_setting_pkey: PRIMARY KEY (id)

Indexes: builtin_wgs_pipeline_setting_pipeline_key_a315120f_like (btree: pipeline_key varchar_pattern_ops)

5.21 checkm2

18 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Checkm2_Accession	character varying(255)	YES			
Name	character varying(255)	YES			
Completeness	character varying(255)	YES			
Contamination	character varying(255)	YES			
Completeness_Model_Used	character varying(255)	YES			
Translation_Table_Used	character varying(255)	YES			
Coding_Density	character varying(255)	YES			
Contig_N50	character varying(255)	YES			
Average_Gene_Length	character varying(255)	YES			
Genome_Size	character varying(255)	YES			
GC_Content	character varying(255)	YES			
Total_Coding_Sequences	character varying(255)	YES			
Total_Contigs	character varying(255)	YES			
Max_Contig_Length	character varying(255)	YES			
Additional_Notes	character varying(255)	YES			
Date_uploaded_c	date	NO			
checkm2_project_id	bigint	YES	FK		WGS_Project.id

Constraints: checkm2_pkey: PRIMARY KEY (id) | checkm2_checkm2_project_id_b6750be1_fk_WGS_Project_id: FOREIGN KEY (checkm2_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: checkm2_Checkm2_Accession_40129c23 (btree: "Checkm2_Accession") | checkm2_Checkm2_Accession_40129c23_like (btree: "Checkm2_Accession" varchar_pattern_ops) | checkm2_checkm2_project_id_b6750be1 (btree: checkm2_project_id)

5.22 Checkm2Upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Checkm2file	character varying(100)	YES			

Constraints: Checkm2Upload_pkey: PRIMARY KEY (id)

5.23 Classification_Table

13 columns; 3 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Class_AccessionNo	character varying(200)	YES			
Class_Chk_Emerging	boolean	NO			
Class_Chk_Satscan	boolean	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
Class_Chk_Serotyping	boolean	NO			
Class_Chk_GHRU_all	boolean	NO			
Class_Chk_GHRU_Neo	boolean	NO			
Class_Chk_Tricycle	boolean	NO			
Class_idNumReferred_id	bigint	NO	FK, UQ		Final_Data.id
Class_Chk_EGASP	boolean	NO			
Class_Chk_Pulsenet	boolean	NO			
Class_Chk_Structured	boolean	NO			
Class_Chk_Tulip	boolean	NO			

Constraints: Classification_Table_pkey: PRIMARY KEY (id) | unique_classification_per_isolate: UNIQUE ("Class_idNumReferred_id") | Classification_Table_Class_idNumReferred__3f17e0d1_fk_Final_Dat: FOREIGN KEY ("Class_idNumReferred_id") REFERENCES public."Final_Data"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: Classification_Table_Class_idNumReferred_id_3f17e0d1 (btree: "Class_idNumReferred_id")

5.24 custom_wgs_pipeline

17 columns; 4 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
name	character varying(150)	NO	UQ		
slug	character varying(170)	NO	UQ		
description	text	NO			
sequencing_type	character varying(20)	NO			
platform	character varying(100)	NO			
category	jsonb	NO			
sheet_name	character varying(120)	NO			
accession_column	character varying(120)	NO			
sample_name_column	character varying(120)	NO			
date_column	character varying(120)	NO			
show_in_upload_center	boolean	NO			
show_in_overview	boolean	NO			
is_active	boolean	NO			
created_at	timestamp with time zone	NO			
updated_at	timestamp with time zone	NO			
created_by_id	integer	YES	FK		auth_user.id

Constraints: custom_wgs_pipeline_name_key: UNIQUE (name) | custom_wgs_pipeline_pkey: PRIMARY KEY (id) | custom_wgs_pipeline_slug_key: UNIQUE (slug) | custom_wgs_pipeline_created_by_id_af63e8f2_fk_auth_user_id: FOREIGN KEY (created_by_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: custom_wgs_pipeline_created_by_id_af63e8f2 (btree: created_by_id) | custom_wgs_pipeline_name_a3da965e_like (btree: name varchar_pattern_ops) | custom_wgs_pipeline_slug_8c3f5733_like (btree: slug varchar_pattern_ops)

5.25 custom_wgs_pipeline_field

13 columns; 4 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
field_key	character varying(120)	NO			
display_label	character varying(150)	NO			
source_column	character varying(150)	NO			
column_aliases	text	NO			
data_type	character varying(20)	NO			
required	boolean	NO			
default_value	character varying(255)	NO			
show_in_table	boolean	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
show_in_detail	boolean	NO			
show_in_export	boolean	NO			
sort_order	integer	NO			
pipeline_id	bigint	NO	FK		custom_wgs_pipeline.id

Constraints: custom_wgs_pipeline_field_sort_order_check: CHECK ((sort_order >= 0)) | custom_wgs_pipeline_field_pipeline_id_field_key_cbe7b481_uniq: UNIQUE (pipeline_id, field_key) | custom_wgs_pipeline_field_pkey: PRIMARY KEY (id) | custom_wgs_pipeline__pipeline_id_360558a9_fk_custom_wg: FOREIGN KEY (pipeline_id) REFERENCES public.custom_wgs_pipeline(id) DEFERRABLE INITIALLY DEFERRED

Indexes: custom_wgs_pipeline_field_field_key_ca715f8a (btree: field_key) | custom_wgs_pipeline_field_field_key_ca715f8a_like (btree: field_key varchar_pattern_ops) | custom_wgs_pipeline_field_pipeline_id_360558a9 (btree: pipeline_id)

5.26 custom_wgs_pipeline_record

14 columns; 6 named constraints; 13 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
accession	character varying(255)	NO			
sample_name	character varying(255)	NO			
match_status	character varying(20)	NO			
match_source	character varying(20)	NO			
values_json	jsonb	NO			
raw_row_json	jsonb	NO			
uploaded_at	timestamp with time zone	NO			
updated_at	timestamp with time zone	NO			
matched_final_data_id	character varying(255)	YES	FK		Final_Data.f_AccessionNo
matched_raw_data_id	character varying(255)	YES	FK		Referred_Data.AccessionNo
pipeline_id	bigint	NO	FK		custom_wgs_pipeline.id
uploaded_by_id	integer	YES	FK		auth_user.id
upload_batch_id	bigint	YES	FK		custom_wgs_pipeline_upload_batch.id

Constraints: custom_wgs_pipeline_record_pkey: PRIMARY KEY (id) | custom_wgs_pipeline__matched_final_data_i_13244890_fk_Final_Dat: FOREIGN KEY (matched_final_data_id) REFERENCES public."Final_Data"("f_AccessionNo") DEFERRABLE INITIALLY DEFERRED | custom_wgs_pipeline__matched_raw_data_id_254a30fb_fk_Referred_: FOREIGN KEY (matched_raw_data_id) REFERENCES public."Referred_Data"("AccessionNo") DEFERRABLE INITIALLY DEFERRED | custom_wgs_pipeline__pipeline_id_b1ddbc96_fk_custom_wg: FOREIGN KEY (pipeline_id) REFERENCES public.custom_wgs_pipeline(id) DEFERRABLE INITIALLY DEFERRED | custom_wgs_pipeline__upload_batch_id_b8262770_fk_custom_wg: FOREIGN KEY (upload_batch_id) REFERENCES public.custom_wgs_pipeline_upload_batch(id) DEFERRABLE INITIALLY DEFERRED | custom_wgs_pipeline__uploaded_by_id_921db27c_fk_auth_user: FOREIGN KEY (uploaded_by_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: custom_wgs__pipelin_a90caf_idx (btree: pipeline_id, accession) | custom_wgs__pipelin_ff5fe9_idx (btree: pipeline_id, match_status) | custom_wgs_pipeline_record_accession_bc272121 (btree: accession) | custom_wgs_pipeline_record_accession_bc272121_like (btree: accession varchar_pattern_ops) | custom_wgs_pipeline_record_match_status_073ff026 (btree: match_status) | custom_wgs_pipeline_record_match_status_073ff026_like (btree: match_status varchar_pattern_ops) | custom_wgs_pipeline_record_matched_final_data_id_13244890 (btree: matched_final_data_id) | custom_wgs_pipeline_record_matched_final_data_id_13244890_like (btree: matched_final_data_id varchar_pattern_ops) | custom_wgs_pipeline_record_matched_raw_data_id_254a30fb (btree: matched_raw_data_id) | custom_wgs_pipeline_record_matched_raw_data_id_254a30fb_like (btree: matched_raw_data_id varchar_pattern_ops) | custom_wgs_pipeline_record_pipeline_id_b1ddbc96 (btree: pipeline_id) | custom_wgs_pipeline_record_upload_batch_id_b8262770 (btree: upload_batch_id) | custom_wgs_pipeline_record_uploaded_by_id_921db27c (btree: uploaded_by_id)

5.27 custom_wgs_pipeline_upload_batch

14 columns; 9 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
file_name	character varying(255)	NO			
sheet_name	character varying(120)	NO			
row_count	integer	NO			
created_count	integer	NO			
updated_count	integer	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
skipped_count	integer	NO			
matched_count	integer	NO			
unmatched_count	integer	NO			
status	character varying(30)	NO			
error_log	text	NO			
uploaded_at	timestamp with time zone	NO			
pipeline_id	bigint	NO	FK		custom_wgs_pipeline.id
uploaded_by_id	integer	YES	FK		auth_user.id

Constraints: custom_wgs_pipeline_upload_batch_created_count_check: CHECK ((created_count >= 0)) | custom_wgs_pipeline_upload_batch_matched_count_check: CHECK ((matched_count >= 0)) | custom_wgs_pipeline_upload_batch_row_count_check: CHECK ((row_count >= 0)) | custom_wgs_pipeline_upload_batch_skipped_count_check: CHECK ((skipped_count >= 0)) | custom_wgs_pipeline_upload_batch_unmatched_count_check: CHECK ((unmatched_count >= 0)) | custom_wgs_pipeline_upload_batch_updated_count_check: CHECK ((updated_count >= 0)) | custom_wgs_pipeline_upload_batch_pkey: PRIMARY KEY (id) | custom_wgs_pipeline__pipeline_id_ed554286_fk_custom_wg: FOREIGN KEY (pipeline_id) REFERENCES public.custom_wgs_pipeline(id) DEFERRABLE INITIALLY DEFERRED | custom_wgs_pipeline__uploaded_by_id_c4d361c7_fk_auth_user: FOREIGN KEY (uploaded_by_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: custom_wgs_pipeline_upload_batch_pipeline_id_ed554286 (btree: pipeline_id) | custom_wgs_pipeline_upload_batch_uploaded_by_id_c4d361c7 (btree: uploaded_by_id)

5.28 DemogsData_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
DemogsDataFile	character varying(100)	YES			

Constraints: DemogsData_upload_pkey: PRIMARY KEY (id)

5.29 django_admin_log

8 columns; 4 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	integer	NO	PK	YES	
action_time	timestamp with time zone	NO			
object_id	text	YES			
object_repr	character varying(200)	NO			
action_flag	smallint	NO			
change_message	text	NO			
content_type_id	integer	YES	FK		django_content_type.id
user_id	integer	NO	FK		auth_user.id

Constraints: django_admin_log_action_flag_check: CHECK ((action_flag >= 0)) | django_admin_log_pkey: PRIMARY KEY (id) | django_admin_log_content_type_id_c4bce8eb_fk_django_co: FOREIGN KEY (content_type_id) REFERENCES public.django_content_type(id) DEFERRABLE INITIALLY DEFERRED | django_admin_log_user_id_c564eba6_fk_auth_user_id: FOREIGN KEY (user_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: django_admin_log_content_type_id_c4bce8eb (btree: content_type_id) | django_admin_log_user_id_c564eba6 (btree: user_id)

5.30 django_content_type

3 columns; 2 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	integer	NO	PK	YES	
app_label	character varying(100)	NO			
model	character varying(100)	NO			

Constraints: django_content_type_app_label_model_76bd3d3b_uniq: UNIQUE (app_label, model) | django_content_type_pkey: PRIMARY KEY (id)

5.31 django_migrations

4 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
app	character varying(255)	NO			
name	character varying(255)	NO			
applied	timestamp with time zone	NO			

Constraints: django_migrations_pkey: PRIMARY KEY (id)

5.32 django_session

3 columns; 1 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
session_key	character varying(40)	NO	PK		
session_data	text	NO			
expire_date	timestamp with time zone	NO			

Constraints: django_session_pkey: PRIMARY KEY (session_key)

Indexes: django_session_expire_date_a5c62663 (btree: expire_date) | django_session_session_key_c0390e0f_like (btree: session_key varchar_pattern_ops)

5.33 Emerging_Filter_Age

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Eme_Age	integer	YES			

Constraints: Emerging_Filter_Age_pkey: PRIMARY KEY (id)

5.34 Emerging_Table

21 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
eme_Spec_Flag	boolean	NO			
eme_Site_Code	character varying(200)	YES			
eme_Accession	character varying(200)	YES			
eme_ReferralDate	character varying(200)	YES			
eme_DateAdmis	character varying(200)	YES			
eme_Diagnosis	character varying(200)	YES			
eme_Diag_ICD	character varying(200)	YES			
eme_ars_Org	character varying(200)	YES			
eme_org_Flag	boolean	NO			
eme_abx_Flag	boolean	NO			
eme_abx_code_pheno	character varying(200)	YES			
eme_abx_Phenotype	character varying(200)	YES			
eme_ars_Pre	character varying(200)	YES			
eme_ars_Post	character varying(200)	YES			
eme_org_Grp	character varying(200)	YES			
eme_org_Genus	character varying(200)	YES			
eme_spec_Num	character varying(200)	YES			
eme_spec_Date	character varying(200)	YES			
eme_spec_Type	character varying(200)	YES			
eme_primary_key_id	bigint	YES	FK		Final_Data.id

Constraints: Emerging_Table_pkey: PRIMARY KEY (id) | Emerging_Table_eme_primary_key_id_e0f72bab_fk_Final_Data_id: FOREIGN KEY (eme_primary_key_id) REFERENCES public."Final_Data"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: Emerging_Table_eme_primary_key_id_e0f72bab (btree: eme_primary_key_id)

5.35 Emerging_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_uploadEme	character varying(100)	YES			

Constraints: Emerging_upload_pkey: PRIMARY KEY (id)

5.36 Final_Data

94 columns; 6 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
f_bat_seq	integer	YES			
f_Hide	boolean	NO			
f_Batch_Code	character varying(255)	NO			
f_Batch_Name	character varying(255)	NO			
f_RefNo	character varying(20)	NO			
f_BatchNo	character varying(255)	NO			
f_Total_batch	character varying(100)	NO			
f_AccessionNo	character varying(255)	NO	UQ		
f_AccessionNoGen	character varying(100)	NO			
f_Date_of_Entry	timestamp with time zone	NO			
f_Date_Modified	timestamp with time zone	NO			
f_SiteCode	character varying(255)	NO			
f_Site_Name	character varying(255)	NO			
f_Referral_Date	date	YES			
f_Patient_ID	character varying(255)	NO			
f_First_Name	character varying(255)	NO			
f_Mid_Name	character varying(255)	NO			
f_Last_Name	character varying(255)	NO			
f_Date_Birth	date	YES			
f_Age	integer	YES			
f_Emerging_Flag_Age	boolean	NO			
f_Sex	character varying(255)	NO			
f_Date_Admis	date	YES			
f_Nosocomial	character varying(255)	NO			
f_Diagnosis	character varying(255)	NO			
f_Diagnosis_ICD10	character varying(255)	NO			
f_Ward	character varying(255)	NO			
f_Ward_Type	character varying(255)	NO			
f_Service_Type	character varying(255)	NO			
f_Spec_Num	character varying(255)	NO			
f_Spec_Date	date	YES			
f_Reason	text	NO			
f_Spec_Emerging	boolean	NO			
f_Growth	character varying(255)	NO			
f_Urine_ColCt	character varying(255)	NO			
f_ampC	character varying(255)	NO			
f_ESBL	character varying(255)	NO			
f_CARB	character varying(255)	NO			
f_MBL	character varying(255)	NO			
f_BL	character varying(255)	NO			
f_MR	character varying(255)	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
f_mecA	character varying(255)	NO			
f_ICR	character varying(255)	NO			
f_OtherResMech	character varying(255)	NO			
f_Site_Pre	character varying(255)	YES			
f_Site_Org	character varying(255)	NO			
f_Site_OrgName	character varying(255)	NO			
f_Site_Pos	character varying(255)	YES			
f_Comments	text	NO			
f_ars_ampC	character varying(255)	NO			
f_ars_ESBL	character varying(255)	NO			
f_ars_CARB	character varying(255)	NO			
f_ars_ECIM	character varying(255)	NO			
f_ars_MCIM	character varying(255)	NO			
f_ars_EC_MCIM	character varying(255)	NO			
f_ars_MBL	character varying(255)	NO			
f_ars_BL	character varying(255)	NO			
f_ars_MR	character varying(255)	NO			
f_ars_mecA	character varying(255)	NO			
f_ars_ICR	character varying(255)	NO			
f_ars_Pre	character varying(255)	YES			
f_ars_Post	character varying(255)	YES			
f_ars_OrgCode	character varying(255)	NO			
f_ars_OrgName	character varying(255)	NO			
f_ars_ct_ctl	character varying(255)	NO			
f_ars_tz_tzl	character varying(255)	NO			
f_ars_cn_cni	character varying(255)	NO			
f_ars_ip_ipi	character varying(255)	NO			
f_ars_reco_Code	character varying(255)	YES			
f_ars_description	text	YES			
f_ars_reco	text	YES			
f_arsp_Encoder	character varying(255)	NO			
f_arsp_Enc_Lic	character varying(100)	NO			
f_arsp_Checker	character varying(255)	NO			
f_arsp_Chec_Lic	character varying(100)	NO			
f_arsp_Verifier	character varying(255)	NO			
f_arsp_Ver_Lic	character varying(100)	NO			
f_arsp_LabManager	character varying(255)	NO			
f_arsp_Lab_Lic	character varying(100)	NO			
f_arsp_Head	character varying(255)	NO			
f_arsp_Head_Lic	character varying(100)	NO			
f_Date_Accomplished_ARSP	date	YES			
f_x_mrse	character varying(255)	NO			
f_x_mrsamrse	character varying(255)	NO			
f_x_entbac	character varying(255)	NO			
f_edta	character varying(255)	NO			
f_Batch_id_id	bigint	YES	FK		Batch_Table.id
f_Spec_Type_id	bigint	YES	FK		SpecimenTypeTable.id
f_Age_Display	character varying(20)	NO			
f_Site_Pos_ed	text	NO			
f_Site_Pre_ed	text	NO			
f_ars_Post_ed	text	NO			
f_ars_Pre_ed	text	NO			

Constraints: Final_Data_f_bat_seq_check: CHECK ((f_bat_seq >= 0)) | Final_Data_f_AccessionNo_key: UNIQUE ("f_AccessionNo") | Final_Data_pkey: PRIMARY KEY (id) | unique_final_accession_batch: UNIQUE ("f_AccessionNo", "f_Batch_Code") | Final_Data_f_Batch_id_id_e5f4c5cf_fk_Batch_Table_id: FOREIGN KEY ("f_Batch_id_id") REFERENCES public."Batch_Table"(id) DEFERRABLE INITIALLY DEFERRED |

Final_Data_f_Spec_Type_id_aaf41bd6_fk_SpecimenTypeTable_id: FOREIGN KEY ("f_Spec_Type_id") REFERENCES public."SpecimenTypeTable"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: Final_Data_f_AccessionNo_e9c0a026_like (btree: "f_AccessionNo" varchar_pattern_ops) | Final_Data_f_Batch_id_id_e5f4c5cf (btree: "f_Batch_id_id") | Final_Data_f_Spec_Type_id_aaf41bd6 (btree: "f_Spec_Type_id")

5.37 FinalAntibiotic_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
FinalAntibioticFile	character varying(100)	YES			

Constraints: FinalAntibiotic_upload_pkey: PRIMARY KEY (id)

5.38 FinalAntibioticEntry

47 columns; 3 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
ab_Site_Org	character varying(100)	YES			
ab_Ret_Org	character varying(100)	YES			
ab_Org_Flag	boolean	NO			
ab_Abx_Flag	boolean	NO			
ab_Abx_Phenotype	character varying(20)	NO			
ab_Disk_Abx	boolean	NO			
ab_Ret_Disk_Abx	boolean	NO			
ab_AccessionNo	character varying(30)	YES			
ab_RefNo	character varying(100)	YES			
ab_Antibiotic	character varying(255)	YES			
ab_Abx_code	character varying(25)	YES			
ab_Abx	character varying(100)	YES			
ab_Disk_value	smallint	YES			
ab_Disk_RIS	character varying(4)	NO			
ab_Disk_enRIS	character varying(4)	NO			
ab_MIC_operand	character varying(4)	YES			
ab_MIC_value	numeric(7,3)	YES			
ab_MIC_RIS	character varying(4)	NO			
ab_MIC_enRIS	character varying(4)	NO			
ab_AlertMIC	boolean	NO			
ab_Alert_val	character varying(30)	YES			
ab_R_breakpoint	character varying(10)	YES			
ab_I_breakpoint	character varying(10)	YES			
ab_SDD_breakpoint	character varying(10)	YES			
ab_S_breakpoint	character varying(10)	YES			
ab_Retest_Antibiotic	character varying(255)	YES			
ab_Retest_Abx_code	character varying(100)	YES			
ab_Retest_Abx	character varying(100)	YES			
ab_Retest_DiskValue	integer	YES			
ab_Retest_Disk_RIS	character varying(4)	NO			
ab_Retest_Disk_enRIS	character varying(4)	NO			
ab_Retest_MIC_operand	character varying(4)	YES			
ab_Retest_MICValue	numeric(7,3)	YES			
ab_Retest_MIC_RIS	character varying(4)	NO			
ab_Retest_MIC_enRIS	character varying(4)	NO			
ab_Retest_AlertMIC	boolean	NO			
ab_Retest_Alert_val	character varying(30)	YES			
ab_Ret_R_breakpoint	character varying(10)	YES			

Column	PostgreSQL type	Null	Key	Identity	Reference
ab_Ret_I_breakpoint	character varying(10)	YES			
ab_Ret_SDD_breakpoint	character varying(10)	YES			
ab_Ret_S_breakpoint	character varying(10)	YES			
ab_Date_uploaded_fd	date	NO			
ab_MICJoined	character varying(7)	YES			
ab_idNum_f_referred_id	bigint	NO	FK		Final_Data.id
ab_Abx_Phenotype_Other	character varying(20)	NO			
ab_eme_specimen	character varying(100)	YES			

Constraints: FinalAntibioticEntry_ab_Disk_value_check: CHECK (("ab_Disk_value" >= 0)) | FinalAntibioticEntry_pkey: PRIMARY KEY (id) | FinalAntibioticEntry_ab_idNum_f_referred__23352223_fk_Final_Dat: FOREIGN KEY ("ab_idNum_f_referred_id") REFERENCES public."Final_Data"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: FinalAntibioticEntry_ab_idNum_f_referred_id_23352223 (btree: "ab_idNum_f_referred_id")

5.39 FinalAntibioticEntry_ab_breakpoints_id

3 columns; 4 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
final_antibioticentry_id	bigint	NO	FK		FinalAntibioticEntry.id
breakpointstable_id	bigint	NO	FK		BreakpointsTable.id

Constraints: FinalAntibioticEntry_ab__final_antibioticentry_id_3968aee4_uniq: UNIQUE (final_antibioticentry_id, breakpointstable_id) | FinalAntibioticEntry_ab_breakpoints_id_pkey: PRIMARY KEY (id) | FinalAntibioticEntry_breakpointstable_id_b2298850_fk_Breakpoint: FOREIGN KEY (breakpointstable_id) REFERENCES public."BreakpointsTable"(id) DEFERRABLE INITIALLY DEFERRED | FinalAntibioticEntry_final_antibioticentr_541d5d3e_fk_FinalAnti: FOREIGN KEY (final_antibioticentry_id) REFERENCES public."FinalAntibioticEntry"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: FinalAntibioticEntry_ab_br_breakpointstable_id_b2298850 (btree: breakpointstable_id) | FinalAntibioticEntry_ab_br_final_antibioticentry_id_541d5d3e (btree: final_antibioticentry_id)

5.40 gambit

15 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Gambit_Accession	character varying(255)	YES			
sample	character varying(255)	YES			
predicted_name	character varying(255)	YES			
predicted_rank	character varying(255)	YES			
predicted_ncbi_id	character varying(255)	YES			
predicted_threshold	character varying(255)	YES			
closest_distance	character varying(255)	YES			
closest_description	character varying(255)	YES			
next_name	character varying(255)	YES			
next_rank	character varying(255)	YES			
next_ncbi_id	character varying(255)	YES			
next_threshold	character varying(255)	YES			
Date_uploaded_g	date	NO			
gambit_project_id	bigint	YES	FK		WGS_Project.id

Constraints: gambit_pkey: PRIMARY KEY (id) | gambit_gambit_project_id_60b79a1a_fk_WGS_Project_id: FOREIGN KEY (gambit_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: gambit_Gambit_Accession_b65a99fb (btree: "Gambit_Accession") | gambit_Gambit_Accession_b65a99fb_like (btree: "Gambit_Accession" varchar_pattern_ops) | gambit_gambit_project_id_60b79a1a (btree: gambit_project_id)

5.41 GambitUpload

2 columns; 1 named constraints; 0 indexes.

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
GambitFile	character varying(100)	YES			

Constraints: GambitUpload_pkey: PRIMARY KEY (id)

5.42 gtdbtk

24 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
GtdbTk_Accession	character varying(255)	YES			
user_genome	character varying(255)	YES			
classification	text	YES			
closest_genome_reference	character varying(255)	YES			
closest_genome_reference_radius	integer	YES			
closest_genome_taxonomy	text	YES			
closest_genome_ani	double precision	YES			
closest_genome_af	double precision	YES			
closest_placement_reference	character varying(255)	YES			
closest_placement_radius	double precision	YES			
closest_placement_taxonomy	text	YES			
closest_placement_ani	double precision	YES			
closest_placement_af	double precision	YES			
pplacer_taxonomy	text	YES			
classification_method	character varying(255)	YES			
note	text	YES			
other_related_references	text	YES			
msa_percent	double precision	YES			
translation_table	integer	YES			
red_value	double precision	YES			
warnings	text	YES			
Date_uploaded_gt	date	NO			
gtdbtk_project_id	bigint	YES	FK		WGS_Project.id

Constraints: gtdbtk_pkey: PRIMARY KEY (id) | gtdbtk_gtdbtk_project_id_76e0ec8e_fk_WGS_Project_id: FOREIGN KEY (gtdbtk_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: gtdbtk_GtdbTk_Accession_b3338e17 (btree: "GtdbTk_Accession") | gtdbtk_GtdbTk_Accession_b3338e17_like (btree: "GtdbTk_Accession" varchar_pattern_ops) | gtdbtk_gtdbtk_project_id_76e0ec8e (btree: gtdbtk_project_id)

5.43 GtdbTkUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
GtdbTkFile	character varying(100)	YES			

Constraints: GtdbTkUpload_pkey: PRIMARY KEY (id)

5.44 home_arsstaff_details

11 columns; 3 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Staff_Name	character varying(100)	YES			
Staff_Designation	character varying(100)	YES			
Staff_Telnum	character varying(128)	YES			
Staff_EmailAdd	character varying(100)	YES			

Column	PostgreSQL type	Null	Key	Identity	Reference
Staff_License	character varying(100)	YES			
Staff_Role	character varying(150)	NO			
User_Account_id	integer	YES	FK, UQ		auth_user.id
Staff_Credentials	character varying(100)	NO			
Is_Default_Head	boolean	NO			
Is_Default_Lab_Manager	boolean	NO			

Constraints: home_arsstaff_details_User_Account_id_key: UNIQUE ("User_Account_id") | home_arsstaff_details_pkey: PRIMARY KEY (id) | home_arsstaff_details_User_Account_id_f48243d2_fk_auth_user_id: FOREIGN KEY ("User_Account_id") REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

5.45 home_fieldmapping

7 columns; 3 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
raw_field	character varying(255)	NO			
mapped_field	character varying(255)	YES			
created_at	timestamp with time zone	NO			
updated_at	timestamp with time zone	NO			
user_id	integer	NO	FK		auth_user.id
is_retest	boolean	NO			

Constraints: home_fieldmapping_pkey: PRIMARY KEY (id) | home_fieldmapping_user_id_raw_field_e6342a63_uniq: UNIQUE (user_id, raw_field) | home_fieldmapping_user_id_6111f8c2_fk_auth_user_id: FOREIGN KEY (user_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

Indexes: home_fieldmapping_user_id_6111f8c2 (btree: user_id)

5.46 home_final_concordancedetail

12 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
accession_no	character varying(100)	NO			
isolate_id	integer	NO			
organism	character varying(255)	NO			
antibiotic	character varying(255)	NO			
site_ris	character varying(3)	NO			
ars_ris	character varying(3)	NO			
report_id	bigint	NO	FK		home_final_concordancereport.id
deviation_code	character varying(1)	YES			
genus_con	character varying(1)	YES			
is_discordant	boolean	NO			
species_con	character varying(1)	YES			

Constraints: home_final_concordancedetail_pkey: PRIMARY KEY (id) | home_final_concordancereport_id_428d5fa5_fk_home_fina: FOREIGN KEY (report_id) REFERENCES public.home_final_concordancereport(id) DEFERRABLE INITIALLY DEFERRED

Indexes: home_final_concordancedetail_report_id_428d5fa5 (btree: report_id)

5.47 home_final_concordancereport

20 columns; 4 named constraints; 6 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
created_at	timestamp with time zone	NO			
total_isolates	integer	NO			
total_pairs	integer	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
concordant_pairs	integer	NO			
vmd	integer	NO			
md	integer	NO			
minor	integer	NO			
ast_concordance_rate	double precision	NO			
batch_id	bigint	YES	FK		Batch_Table.id
created_by_id	integer	YES	FK		auth_user.id
critical_deviation	integer	NO			
critical_deviation_rate	double precision	NO			
total_deviation	integer	NO			
total_deviation_rate	double precision	NO			
final_data_id	bigint	YES	FK		Final_Data.id
genus_match	integer	NO			
genus_rate	double precision	NO			
species_match	integer	NO			
species_rate	double precision	NO			

Constraints: home_final_concordancereport_pkey: PRIMARY KEY (id) | home_final_concordancereport_batch_id_d6ffbf20_fk_Batch_Tab: FOREIGN KEY (batch_id) REFERENCES public."Batch_Table"(id) DEFERRABLE INITIALLY DEFERRED | home_final_concordancereport_created_by_id_d9a0cf3a_fk_auth_user: FOREIGN KEY (created_by_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED | home_final_concordancereport_final_data_id_466017ec_fk_Final_Data: FOREIGN KEY (final_data_id) REFERENCES public."Final_Data"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: concordance_created_idx (btree: created_at DESC) | home_final_concordancereport_batch_id_d6ffbf20 (btree: batch_id) | home_final_concordancereport_created_by_id_d9a0cf3a (btree: created_by_id) | home_final_concordancereport_final_data_id_466017ec (btree: final_data_id) | unique_accession_concordance (UNIQUE btree: final_data_id) WHERE (final_data_id IS NOT NULL) | unique_batch_concordance (UNIQUE btree: batch_id) WHERE (batch_id IS NOT NULL)

5.48 home_nonworkingday

5 columns; 2 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
date	date	NO	UQ		
description	character varying(255)	NO			
applies_to	character varying(10)	NO			
created_at	timestamp with time zone	NO			

Constraints: home_nonworkingday_date_key: UNIQUE (date) | home_nonworkingday_pkey: PRIMARY KEY (id)

5.49 home_organism_list

16 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Whonet_Org_Code	character varying(20)	NO	UQ		
Replaced_by	character varying(20)	YES			
Organism	character varying(255)	NO			
Organism_Type	character varying(5)	YES			
Family_Code	character varying(20)	YES			
Genus_Group	character varying(50)	YES			
Genus_Code	character varying(20)	YES			
Species_Group	character varying(50)	YES			
Serovar_Group	character varying(50)	YES			
Kingdom	character varying(100)	YES			
Phylum	character varying(100)	YES			
Class	character varying(100)	YES			
Order	character varying(100)	YES			
Family	character varying(100)	YES			

Column	PostgreSQL type	Null	Key	Identity	Reference
Genus	character varying(100)	YES			

Constraints: home_organism_list_Whonet_Org_Code_key: UNIQUE ("Whonet_Org_Code") | home_organism_list_pkey: PRIMARY KEY (id)

Indexes: home_organism_list_Whonet_Org_Code_cf78fc7f_like (btree: "Whonet_Org_Code" varchar_pattern_ops)

5.50 home_pheno_upload_pre

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_Pheno_pre	character varying(100)	YES			

Constraints: home_pheno_upload_pre_pkey: PRIMARY KEY (id)

5.51 home_reco_item_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_reco_desc	character varying(100)	YES			

Constraints: home_reco_item_upload_pkey: PRIMARY KEY (id)

5.52 home_recommendation

3 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Reco_Code	character varying(100)	YES			
Reco_Details	text	YES			

Constraints: home_recommendation_pkey: PRIMARY KEY (id)

5.53 home_sitedata

19 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
SiteCode	character varying(3)	NO			
SiteName	character varying(155)	NO			
Site_Address	character varying(255)	NO			
Site_Lab_Head	character varying(255)	NO			
Site_Lab_Head_Contact	character varying(100)	NO			
Site_Lab_Head_Designation	character varying(255)	NO			
Site_Lab_Head_Email	character varying(254)	NO			
Site_MedTech	character varying(255)	NO			
Site_MedTech_Contact	character varying(100)	NO			
Site_MedTech_Designation	character varying(255)	NO			
Site_MedTech_Email	character varying(254)	NO			
Site_Med_Ctr_Chief	character varying(255)	NO			
Site_Med_Ctr_Chief_Contact	character varying(100)	NO			
Site_Med_Ctr_Chief_Designation	character varying(255)	NO			
Site_Med_Ctr_Chief_Email	character varying(254)	NO			
Site_Lab_Head_Credentials	character varying(100)	NO			
Site_MedTech_Credentials	character varying(100)	NO			
Site_Med_Ctr_Chief_Credentials	character varying(100)	NO			

Constraints: home_sitedata_pkey: PRIMARY KEY (id)

5.54 home_tatstep

11 columns; 5 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
step_type	character varying(500)	NO			
step_owner	character varying(10)	NO			
step_days_count	integer	YES			
within_tat	boolean	YES			
date_received	date	YES			
date_finished	date	YES			
remarks	text	YES			
performed_by_id	bigint	YES	FK		home_arsstaff_details.id
tat_id	bigint	NO	FK		TATform.id
step_config_id	bigint	NO	FK		home_tatstepconfig.id

Constraints: home_tatstep_step_days_count_check: CHECK ((step_days_count >= 0)) | home_tatstep_pkey: PRIMARY KEY (id) | home_tatstep_performed_by_id_598e377c_fk_home_arss: FOREIGN KEY (performed_by_id) REFERENCES public.home_arsstaff_details(id) DEFERRABLE INITIALLY DEFERRED | home_tatstep_step_config_id_6d1898ea_fk_home_tatstepconfig_id: FOREIGN KEY (step_config_id) REFERENCES public.home_tatstepconfig(id) DEFERRABLE INITIALLY DEFERRED | home_tatstep_tat_id_2d827def_fk_TATform_id: FOREIGN KEY (tat_id) REFERENCES public."TATform"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: home_tatstep_performed_by_id_598e377c (btree: performed_by_id) | home_tatstep_step_config_id_6d1898ea (btree: step_config_id) | home_tatstep_tat_id_2d827def (btree: tat_id)

5.55 home_tatstepconfig

5 columns; 3 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
step_type	character varying(500)	NO			
step_owner	character varying(255)	NO			
target_days	integer	NO			
order	integer	NO			

Constraints: home_tatstepconfig_order_check: CHECK (("order" >= 0)) | home_tatstepconfig_target_days_check: CHECK ((target_days >= 0)) | home_tatstepconfig_pkey: PRIMARY KEY (id)

5.56 mlst

14 columns; 2 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
MIst_Accession	character varying(255)	YES			
name	character varying(255)	YES			
scheme	character varying(255)	YES			
mlst	character varying(255)	YES			
allele1	character varying(255)	YES			
allele2	character varying(255)	YES			
allele3	character varying(255)	YES			
allele4	character varying(255)	YES			
allele5	character varying(255)	YES			
allele6	character varying(255)	YES			
allele7	character varying(255)	YES			
Date_uploaded_m	date	NO			
mlst_project_id	bigint	YES	FK		WGS_Project.id

Constraints: mlst_pkey: PRIMARY KEY (id) | mlst_mlst_project_id_24fb710b_fk_WGS_Project_id: FOREIGN KEY (mlst_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: mlst_Mlst_Accession_8f6e7032 (btree: "Mlst_Accession") | mlst_Mlst_Accession_8f6e7032_like (btree: "Mlst_Accession" varchar_pattern_ops) | mlst_mlst_project_id_24fb710b (btree: mlst_project_id)

5.57 MlstUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Mlstfile	character varying(100)	YES			

Constraints: MlstUpload_pkey: PRIMARY KEY (id)

5.58 Organism_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_uploadOrg	character varying(100)	YES			

Constraints: Organism_upload_pkey: PRIMARY KEY (id)

5.59 Pheno_upload_Post

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_Pheno_post	character varying(100)	YES			

Constraints: Pheno_upload_Post_pkey: PRIMARY KEY (id)

5.60 Phenotype_Post

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Post_Phenotypes	character varying(255)	YES			

Constraints: Phenotype_Post_pkey: PRIMARY KEY (id)

5.61 Phenotype_Pre

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Pre_Phenotypes	character varying(255)	YES			

Constraints: Phenotype_Pre_pkey: PRIMARY KEY (id)

5.62 RawAntibiotic_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
RawAntibioticFile	character varying(100)	YES			

Constraints: RawAntibiotic_upload_pkey: PRIMARY KEY (id)

5.63 Recommendation_items

3 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	

Column	PostgreSQL type	Null	Key	Identity	Reference
RecoCode	character varying(100)	YES			
Description	text	YES			

Constraints: Recommendation_items_pkey: PRIMARY KEY (id)

5.64 Referred_Data

107 columns; 6 named constraints; 3 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
bat_seq	integer	YES			
Hide	boolean	NO			
Copy_data	boolean	NO			
Batch_Name	character varying(255)	NO			
Batch_Code	character varying(255)	NO			
Date_of_Entry	timestamp with time zone	NO			
Date_Modified	timestamp with time zone	NO			
RefNo	character varying(20)	YES			
BatchNo	character varying(255)	NO			
Total_batch	character varying(100)	NO			
AccessionNo	character varying(255)	NO	UQ		
AccessionNoGen	character varying(100)	NO			
Default_Year	date	YES			
SiteCode	character varying(255)	NO			
Site_Name	character varying(255)	NO			
Referral_Date	date	YES			
Patient_ID	character varying(255)	NO			
First_Name	character varying(255)	NO			
Mid_Name	character varying(255)	NO			
Last_Name	character varying(255)	NO			
Date_Birth	date	YES			
Age	integer	YES			
Sex	character varying(255)	NO			
Date_Admis	date	YES			
Nosocomial	character varying(255)	NO			
Diagnosis	character varying(255)	NO			
Diagnosis_ICD10	character varying(255)	NO			
Ward	character varying(255)	NO			
Ward_Type	character varying(255)	NO			
Service_Type	character varying(255)	NO			
Spec_Num	character varying(255)	NO			
Spec_Date	date	YES			
Reason	text	NO			
Growth	character varying(255)	NO			
Urine_ColCt	character varying(255)	NO			
ampC	character varying(255)	NO			
ESBL	character varying(255)	NO			
CARB	character varying(255)	NO			
MBL	character varying(255)	NO			
BL	character varying(255)	NO			
MR	character varying(255)	NO			
mecA	character varying(255)	NO			
ICR	character varying(255)	NO			
OtherResMech	character varying(255)	NO			
Site_Pre	character varying(255)	NO			
Site_Org	character varying(255)	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
Site_OrgName	character varying(255)	YES			
Site_Pos	character varying(255)	NO			
Comments	text	YES			
ars_ampC	character varying(255)	NO			
ars_ESBL	character varying(255)	NO			
ars_CARB	character varying(255)	NO			
ars_ECIM	character varying(255)	NO			
ars_MCIM	character varying(255)	NO			
ars_EC_MCIM	character varying(255)	NO			
ars_MBL	character varying(255)	NO			
ars_BL	character varying(255)	NO			
ars_MR	character varying(255)	NO			
ars_mecA	character varying(255)	NO			
ars_ICR	character varying(255)	NO			
ars_Pre	character varying(255)	NO			
ars_Post	character varying(255)	NO			
ars_OrgCode	character varying(255)	NO			
ars_OrgName	character varying(255)	NO			
ars_ct_ctl	character varying(255)	NO			
ars_tz_tzl	character varying(255)	NO			
ars_cn_cni	character varying(255)	NO			
ars_ip_ipi	character varying(255)	NO			
ars_reco_Code	character varying(255)	YES			
ars_description	text	YES			
ars_reco	text	YES			
SiteName	character varying(255)	NO			
Status	character varying(100)	NO			
Month_Date	date	YES			
Day_Date	date	YES			
Year_Date	date	YES			
RefDate	date	YES			
Start_AccNo	integer	YES			
End_AccNo	integer	YES			
No_Isolates	integer	YES			
Concordance_Check	character varying(255)	NO			
Concordance_by	character varying(255)	NO			
Concordance_by_Initials	character varying(255)	NO			
abx_code	character varying(25)	NO			
arsp_Encoder	character varying(255)	YES			
arsp_Enc_Lic	character varying(100)	YES			
arsp_Checker	character varying(255)	YES			
arsp_Chec_Lic	character varying(100)	YES			
arsp_Verifier	character varying(255)	YES			
arsp_Ver_Lic	character varying(100)	YES			
arsp_LabManager	character varying(255)	YES			
arsp_Lab_Lic	character varying(100)	YES			
arsp_Head	character varying(255)	YES			
arsp_Head_Lic	character varying(100)	YES			
Date_Accomplished_ARSP	date	YES			
x_mrse	character varying(255)	NO			
x_mrsamrse	character varying(255)	NO			
x_entbac	character varying(255)	NO			
edta	character varying(255)	NO			
Batch_id_id	bigint	YES	FK		Batch_Table.id
Spec_Type_id	bigint	YES	FK		SpecimenTypeTable.id

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
Age_Display	character varying(20)	NO			
Site_Pos_ed	text	NO			
Site_Pre_ed	text	NO			
ars_Post_ed	text	NO			
ars_Pre_ed	text	NO			

Constraints: Referred_Data_bat_seq_check: CHECK ((bat_seq >= 0)) | Referred_Data_AccessionNo_key: UNIQUE ("AccessionNo") | Referred_Data_pkey: PRIMARY KEY (id) | unique_accession_batch: UNIQUE ("AccessionNo", "Batch_Code") | Referred_Data_Batch_id_id_2c0dc7b7_fk_Batch_Table_id: FOREIGN KEY ("Batch_id_id") REFERENCES public."Batch_Table"(id) DEFERRABLE INITIALLY DEFERRED | Referred_Data_Spec_Type_id_35f213b0_fk_SpecimenTypeTable_id: FOREIGN KEY ("Spec_Type_id") REFERENCES public."SpecimenTypeTable"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: Referred_Data_AccessionNo_7f8b786d_like (btree: "AccessionNo" varchar_pattern_ops) | Referred_Data_Batch_id_id_2c0dc7b7 (btree: "Batch_id_id") | Referred_Data_Spec_Type_id_35f213b0 (btree: "Spec_Type_id")

5.65 Referred_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
ReferredDataFile	character varying(100)	YES			

Constraints: Referred_upload_pkey: PRIMARY KEY (id)

5.66 sample_information

25 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
sample_accession	character varying(255)	YES			
batch_code	character varying(255)	YES			
sample_name	character varying(255)	YES			
status	character varying(255)	YES			
emerging	boolean	NO			
structured	boolean	NO			
satscan	boolean	NO			
serotyping	boolean	NO			
ghru	boolean	NO			
egasp	boolean	NO			
tricycle	boolean	NO			
pulsenet	boolean	NO			
tulip	boolean	NO			
Date_uploaded_si	date	NO			
sample_project_id	bigint	YES	FK		WGS_Project.id
ghru_all	boolean	NO			
ghru_neo	boolean	NO			
DNA_extraction	boolean DEFAULT false	NO			
dna_performed_by	character varying(255)	YES			
dna_performed_date	date	YES			
library_performed_by	character varying(255)	YES			
library_performed_date	date	YES			
library_preparation	boolean DEFAULT false	NO			
sequencing_platform	boolean DEFAULT false	NO			

Constraints: sample_information_pkey: PRIMARY KEY (id) | sample_information_sample_project_id_7ad83fbd_fk_WGS_Project_id: FOREIGN KEY (sample_project_id) REFERENCES public."WGS_Project"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: sample_information_sample_project_id_7ad83fbd (btree: sample_project_id)

5.67 SampleInfoUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
sampleinfo	character varying(100)	YES			

Constraints: SampleInfoUpload_pkey: PRIMARY KEY (id)

5.68 SiteCode_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_uploadSite	character varying(100)	YES			

Constraints: SiteCode_upload_pkey: PRIMARY KEY (id)

5.69 Specimen_upload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
File_uploadSpec	character varying(100)	YES			

Constraints: Specimen_upload_pkey: PRIMARY KEY (id)

5.70 SpecimenTypeTable

6 columns; 3 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Emerging_Spec_Flag	boolean	NO			
Specimen_code	character varying(4)	YES	UQ		
Specimen_name	character varying(100)	YES			
Specimen_Grp_Name	character varying(100)	YES			
Specimen_Code_Grp_id	bigint	YES	FK		SpecimenTypeTable.id

Constraints: SpecimenTypeTable_Specimen_code_key: UNIQUE ("Specimen_code") | SpecimenTypeTable_pkey: PRIMARY KEY (id) | SpecimenTypeTable_Specimen_Code_Grp_id_22b1f812_fk_SpecimenT: FOREIGN KEY ("Specimen_Code_Grp_id") REFERENCES public."SpecimenTypeTable"(id) DEFERRABLE INITIALLY DEFERRED

Indexes: SpecimenTypeTable_Specimen_Code_Grp_id_22b1f812 (btree: "Specimen_Code_Grp_id") | SpecimenTypeTable_Specimen_code_14775a5d_like (btree: "Specimen_code" varchar_pattern_ops)

5.71 TATform

19 columns; 7 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
tat_SiteCode	character varying(10)	NO			
tat_Batch_Code	character varying(255)	NO			
tat_Referral_Date	date	YES			
tat_BatchNumber	character varying(255)	NO			
tat_Total_Batch	character varying(100)	NO			
tat_Num_Isolate	integer	YES			
tat_Batch_Location	character varying(20)	NO			
tat_Target_Days	integer	YES			
tat_Date_Released	date	YES			
tat_Running_TAT	integer	YES			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
tat_Final_TAT	integer	YES			
tat_Status_Release	character varying(20)	NO			
tat_Remarks	text	YES			
tat_Date_Last_Update	date	YES			
tat_Batch_Isolates_id	bigint	NO	FK, UQ		Batch_Table.id
tat_Scanning_final	boolean	NO			
tat_Scanning_raw	boolean	NO			
tat_Scanning_ws	boolean	NO			

Constraints: TATform_tat_Final_TAT_check: CHECK (("tat_Final_TAT" >= 0)) | TATform_tat_Num_Isolate_check: CHECK (("tat_Num_Isolate" >= 0)) | TATform_tat_Running_TAT_check: CHECK (("tat_Running_TAT" >= 0)) | TATform_tat_Target_Days_check: CHECK (("tat_Target_Days" >= 0)) | TATform_pkey: PRIMARY KEY (id) | TATform_tat_Batch_Isolates_id_key: UNIQUE ("tat_Batch_Isolates_id") | TATform_tat_Batch_Isolates_id_464d7e8d_fk_Batch_Table_id: FOREIGN KEY ("tat_Batch_Isolates_id") REFERENCES public."Batch_Table"(id) DEFERRABLE INITIALLY DEFERRED

5.72 TATStepConfigUpload

2 columns; 1 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
tat_file	character varying(100)	YES			

Constraints: TATStepConfigUpload_pkey: PRIMARY KEY (id)

5.73 UserProfile

3 columns; 3 named constraints; 0 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
Middle_Name	character varying(150)	NO			
user_id	integer	NO	FK, UQ		auth_user.id

Constraints: UserProfile_pkey: PRIMARY KEY (id) | UserProfile_user_id_key: UNIQUE (user_id) | UserProfile_user_id_d5588aad_fk_auth_user_id: FOREIGN KEY (user_id) REFERENCES public.auth_user(id) DEFERRABLE INITIALLY DEFERRED

5.74 wgs_app_gambitdisplayconfig

3 columns; 2 named constraints; 1 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
field_name	character varying(100)	NO	UQ		
show	boolean	NO			

Constraints: wgs_app_gambitdisplayconfig_field_name_key: UNIQUE (field_name) | wgs_app_gambitdisplayconfig_pkey: PRIMARY KEY (id)

Indexes: wgs_app_gambitdisplayconfig_field_name_d2c32c4c_like (btree: field_name varchar_pattern_ops)

5.75 WGS_Project

18 columns; 2 named constraints; 2 indexes.

Column	PostgreSQL type	Null	Key	Identity	Reference
id	bigint	NO	PK	YES	
WGS_BactScout_Acc	character varying(255)	YES			
WGS_BactScoutSummary	boolean	NO			
WGS_Gambit_Acc	character varying(255)	YES			
WGS_GambitSummary	boolean	NO			
WGS_Mlst_Acc	character varying(255)	YES			
WGS_MlstSummary	boolean	YES			
WGS_Checkm2_Acc	character varying(255)	YES			
WGS_Checkm2Summary	boolean	NO			

REFERRED_DB Database Schema

Column	PostgreSQL type	Null	Key	Identity	Reference
WGS_Assembly_Acc	character varying(255)	YES			
WGS_AssemblySummary	boolean	NO			
WGS_Amrfinder_Acc	character varying(255)	YES			
WGS_AmrfinderSummary	boolean	NO			
Ref_Accession_id	character varying(255)	YES	FK		Final_Data.f_AccessionNo
WGS_GtdbTkSummary	boolean	NO			
WGS_GtdbTk_Acc	character varying(255)	YES			
WGS_SampleInfoSummary	boolean	NO			
WGS_SampleInfo_Acc	character varying(255)	YES			

Constraints: WGS_Project_pkey: PRIMARY KEY (id) | WGS_Project_Ref_Accession_id_0dcf3ffa_fk_Final_Dat: FOREIGN KEY ("Ref_Accession_id") REFERENCES public."Final_Data"("f_AccessionNo") DEFERRABLE INITIALLY DEFERRED

Indexes: WGS_Project_Ref_Accession_id_0dcf3ffa (btree: "Ref_Accession_id") | WGS_Project_Ref_Accession_id_0dcf3ffa_like (btree: "Ref_Accession_id" varchar_pattern_ops)

6. Schema Notes

- The supplied dump is schema-only, so this document does not include row counts or sample records.
- PostgreSQL identity definitions are shown as IDENTITY in the detailed table catalog.
- Many-to-many junction tables are represented as ordinary tables with their composite UNIQUE constraints and foreign keys.
- Application-level semantics that are not encoded in SQL (for example, business descriptions of each field) cannot be determined from the schema dump alone.
- The organizational grouping used in this document is an interpretation based on table names and Django conventions; the exact Django app/model mapping is best confirmed from the project models and migrations.